

CellEnrich Manual

2023.11.11 @hainct

Introduction

Testing environment:

1. Ubuntu 22.04.2 LTS
2. R version 4.3.0
3. RStudio 2023.03.0+386 “Cherry Blossom”
4. Rtools 40

NOTE: on a fresh installation, users may need to install some required interpreter compilers for the system to install other R packages:

- C++ compiler
- gfortran compiler (FYI: tips for [MAC](#) users or [other OS](#))
- Seurat >= version 5.0.0 is REQUIRED

Install Dependent Packages:

- Also, some package needs **install from sources which needs compilation**, and it needs proper version of [Rtools](#) installed
- It's recommended to install [devtools](#), [BiocManager](#) packages first, to install following packages in Github / Bioconductor (not in CRAN).

```
library(devtools)
library(BiocManager)
```

Bioconductor Packages:

1. SingleCellExperiment
2. Scraper

```
BiocManager::install('SingleCellExperiment')
BiocManager::install('scraper')
```

Other Packages:

1. Seurat
2. farver
3. ggbiplot
4. waiter

```
install.packages('Seurat') # version >= 5.0.0 is REQUIRED!
install.packages('farver')
devtools::install_github('JohnCoene/waiter')
devtools::install_github('vqv/ggbiplot')
```

After installing these packages, we can install CellEnrich now.

Running this script will install most of the other dependent packages.

```
install_github('noobCoding/CellEnrich')
library(CellEnrich)
```

Download Data

CellEnrich provides example data and genesets.

Using this function, users can download them in the current working directory.

```
CellEnrich::DownloadData(type='test') # for available human
related data & genesets
```

```
## Use type='all' for all available data
```

Here is a list of data.

- Human Gene-sets
 - Reactome_2022.RData
 - WikiPathways_2021_Human.RData
 - KEGG_2021_Human.RData
 - humanGO.RData
 - humanGOBP.RData
 - humanGOCC.RData
 - humanGOMF.RData
- Mouse Gene-sets
 - WikiPathways_2019_Mouse.RData
 - KEGG_2021_Mouse.RData
 - mouseGO.RData
 - mouseGOBP.RData
 - mouseGOCC.RData
 - mouseGOMF.RData
- Human Data
 - pbmcClustInfo.RData ([Seurat PBMC guide](#))
 - pbmcData.RData
 - [DuoClustering\(link\)](#)
 - koh.RData
 - kohInfo.RData
- Mouse Data
 - [Seurat MCA guide](#)

Tutorial with PBMC 3k data

This tutorial will use PBMC analysis data from [Seurat](#).

To get a complete description of the experiment, refer [here](#)

We assume that all dependent packages are already installed. See **How to install CellEnrich** document.

Required Data

Pre-processed data is included in GitHub [repository](#).

| Filename | Parameter | Type |
|---------------------|-----------|-----------|
| pbmcData.RData | CountData | dgCMatrix |
| pbmcClustInfo.RData | GroupInfo | Character |

Run these codes, users can see an interactive shiny page.

```
# Download data, if not downloaded
CellEnrich::DownloadData()

load("pbmcData.RData")
load("pbmcClustInfo.RData")

CountData <- pbmcData
GroupInfo <- pbmcClustInfo

# This will run CellEnrich
CellEnrich(CountData, GroupInfo)
```

Result

Set options before starting CellEnrich.

Possible options are :

- Strategy for Differential Expressed Genes.
- Plot Options
- Gene-set and their size cutoff to be used in the analysis.

After a few minutes, the analysis result will appear.

Test environment used:

- OS : Ubuntu 22.04 LTS
- CPU: iAMD® Ryzen 9 5950x 16-core processor × 32
- RAM: 64GB

CellEnrich consists of 5 modules.

2. Scatter Area

- The top-left plot is a scatter plot created with PCA, TSNE or UMAP and ggplot2 that the user selected in the option.
- The top-right plot is a histogram plot created with high charter to see a distribution of Group / Cell labels.
- The bottom-left plot shows the significance level (darker is more significant) of enriched pathways dominant in each cell of a group.
- The bottom-right plot shows the significance level of enriched pathways in the whole data.



- Users can use the designed buttons to change the enriched criteria:
 1. ODDS RATIO button will colorize the enriched cell in each group with the **highest** odds ratio.
 2. FREQUENCY button will colorize frequently enriched cells in each group.
- Users also can use the designed buttons to download analyzed figures and results.

3. Pathways Module

This pathway module will show significant pathways for each group based on the Frequency (of pathways in a group), Size, OddsRatio and Q-value (of pathways) in table format.

For each group, **only 1 pathway** can be selected to show in the scatter plot at a time.

B

SELECT

| Geneset | Frequency | Size | Qvalue | OddRatio |
|---|-----------|------|--------|----------|
| ANTIGEN PROCESSING AND PRESENTATION OF EXOGENOUS PEPTIDE ANTIGEN VIA MHC CLASS II | 343 | 27 | 20 | 79.2278 |
| ER TO GOLGI TRANSPORT VESICLE MEMBRANE | 343 | 51 | 20 | 79.0949 |
| ENDOCYTIC VESICLE MEMBRANE | 340 | 45 | 20 | 77.6131 |
| POSITIVE REGULATION OF T CELL ACTIVATION | 336 | 33 | 20 | 70.1417 |
| TRANSPORT VESICLE MEMBRANE | 315 | 27 | 20 | 42.898 |
| ADAPTIVE IMMUNE RESPONSE | 299 | 159 | 20 | 15.4034 |
| CLATHRIN-COATED ENDOCYTIC VESICLE MEMBRANE | 290 | 42 | 20 | 11.781 |
| B CELL RECEPTOR SIGNALING PATHWAY | 98 | 42 | 20 | 9.3323 |
| POSITIVE REGULATION OF T CELL MEDIATED CYTOTOXICITY | 291 | 21 | 20 | 8.6265 |
| B CELL ACTIVATION | 52 | 34 | 20 | 8.2109 |

Previous

1

2

3

Next

CD8 T

SELECT

| Geneset | Frequency | Size | Qvalue | OddRatio |
|---|-----------|------|---------|----------|
| NATURAL KILLER CELL MEDIATED CYTOTOXICITY | 28 | 22 | 13.7615 | 5.4114 |
| T CELL RECEPTOR SIGNALING PATHWAY | 40 | 84 | 15.9546 | 4.6352 |
| IMMUNOLOGICAL SYNAPSE | 98 | 42 | 20 | 3.1694 |
| RRNA PROCESSING | 257 | 146 | 5.8533 | 2.9291 |
| MHC CLASS I PROTEIN BINDING | 33 | 17 | 6.8354 | 2.7643 |
| T CELL ACTIVATION | 65 | 55 | 11.5 | 2.7447 |
| VESICLE | 236 | 129 | 3.9055 | 1.8295 |
| PROTEOLYSIS INVOLVED IN PROTEIN CATABOLIC PROCESS | 101 | 44 | 5.3877 | 1.7347 |
| RIBOSOMAL LARGE SUBUNIT ASSEMBLY | 245 | 17 | 2.0803 | 1.5907 |
| UNFOLDED PROTEIN BINDING | 117 | 94 | 3.5068 | 1.5049 |

Previous

1

2

Next

CD14+ Mono

SELECT

| Geneset | Frequency | Size | Qvalue | OddRatio |
|---|-----------|------|--------|----------|
| SECRETORY GRANULE LUMEN | 477 | 100 | 20 | 35.0065 |
| FICOLIN-1-RICH GRANULE LUMEN | 476 | 118 | 20 | 30.1316 |
| TERTIARY GRANULE LUMEN | 444 | 41 | 20 | 28.9028 |
| ANTIGEN PROCESSING AND PRESENTATION OF EXOGENOUS PEPTIDE ANTIGEN VIA MHC CLASS II | 442 | 27 | 20 | 17.4529 |
| COLLAGEN-CONTAINING EXTRACELLULAR MATRIX | 358 | 134 | 20 | 14.3061 |
| ANTIGEN PROCESSING AND PRESENTATION | 448 | 39 | 20 | 13.0168 |
| POSITIVE REGULATION OF INFLAMMATORY RESPONSE | 253 | 78 | 20 | 8.9201 |
| MHC CLASS II PROTEIN COMPLEX | 402 | 15 | 20 | 7.8069 |
| ENDOCYTIC VESICLE MEMBRANE | 391 | 45 | 20 | 7.2544 |
| NEUTROPHIL CHEMOTAXIS | 164 | 47 | 20 | 6.5747 |

Previous

1

2

3

4

5

...

10

Next

DC

SELECT

| Geneset | Frequency | Size | Qvalue | OddRatio |
|--|-----------|------|---------|----------|
| ER TO GOLGI TRANSPORT VESICLE MEMBRANE | 31 | 51 | 7.7155 | 28.823 |
| COPII-COATED ER TO GOLGI TRANSPORT VESICLE | 5 | 34 | 5.7082 | 25.5263 |
| ENDOCYTIC VESICLE MEMBRANE | 30 | 45 | 10.5129 | 24.7688 |
| PROTON TRANSMEMBRANE TRANSPORT | 31 | 80 | 5.2115 | 18.3829 |
| LYSOSOMAL MEMBRANE | 27 | 294 | 12.7768 | 18.342 |
| TRANSPORT VESICLE MEMBRANE | 26 | 27 | 12.8556 | 17.1138 |
| LYSOSOME | 25 | 369 | 12.414 | 15.4234 |
| CLATHRIN-COATED ENDOCYTIC VESICLE MEMBRANE | 28 | 42 | 10.3106 | 15.3559 |
| OXIDATIVE PHOSPHORYLATION | 31 | 20 | 4.263 | 14.8654 |
| RESPONSE TO ENDOPLASMIC RETICULUM STRESS | 4 | 73 | 3.655 | 14.3516 |

Previous

1

2

3

4

5

...

11

Next

User chosen Pathways

CELLULAR RESPONSE TO TUMOR NECROSIS FACTOR @NK

RRNA PROCESSING @Naive CD4 T

COPII-COATED ER TO GOLGI TRANSPORT VESICLE @DC

RESPIRASOME @Memory CD4 T

TERTIARY GRANULE LUMEN @FCGR3A+ Mono

TERTIARY GRANULE LUMEN @CD14+ Mono

PLOT THE SELECTED PATHWAYS



















CLEAR LIST

PLEASE SWITCH ANY ELEMENT'S POSITION TO ACTIVATE THE PLOT BUTTON

NOTE: please switch any element’s position to activate the plot button

After clicking the ‘PLOT THE SELECTED PATHWAYS’ button, the top cell enriched by the selected pathways will be shown in the Scatter & Bar session.


The corresponding legend will also be displayed. And it can be downloaded separately.

| Scale | | Pathway | Group |
|--|---|---|--------------|
| High  Low | |  B CELL ACTIVATION | B |
| High  Low |  | POSITIVE REGULATION OF INFLAMMATORY RESPONSE | CD14+ Mono |
| High  Low |  | NATURAL KILLER CELL MEDIATED CYTOTOXICITY | CD8 T |
| High  Low |  | COPII-COATED ER TO GOLGI TRANSPORT VESICLE | DC |
| High  Low |  | SPECIFIC GRANULE MEMBRANE | FCGR3A+ Mono |
| High  Low |  | T CELL ACTIVATION | Memory CD4 T |
| High  Low |  | CELLULAR RESPONSE TO INTERLEUKIN-4 | Naive CD4 T |
| High  Low |  | MONOCYTE CHEMOTAXIS | NK |
| High  Low |  | PLATELET FORMATION | Platelet |


4. BiPlot & HeatMap Module

This section shows the BiPlot presentation of the correlation between categories and top significant pathways.

Biplot between top pathways and groups



Users can check the BiPlot under the perspective of **FREQUENCY** or **ODDS RATIO**. This will also update the corresponding HeatMap.



5. Marker Module

MarkerGenes

DE from each Cell specific

| Group | Top | genes | FDR | | | | | |
|----------|-----|--------|-----|---|---|-----|----|------|
| All | All | All | All | | | | | |
| 1 | 1 | CD3D | 0 | | | | | |
| 1 | 1 | RPS27A | 0 | | | | | |
| 1 | 1 | RPS27 | 0 | | | | | |
| 1 | 3 | LTB | 0 | | | | | |
| 1 | 5 | RPL9 | 0 | | | | | |
| 1 | 2 | RPL21 | 0 | | | | | |
| 1 | 3 | RPS3A | 0 | | | | | |
| 1 | 2 | LDHB | 0 | | | | | |
| 1 | 4 | RPS25 | 0 | | | | | |
| 1 | 2 | CD3E | 0 | | | | | |
| Previous | 1 | 2 | 3 | 4 | 5 | ... | 36 | Next |

DE - Pathway from each Cell specific

| genes | Count | Group | Top | FDR | | | | |
|----------|-------|-------|-----|-----|---|-----|----|------|
| All | All | All | All | All | | | | |
| ACTN1 | 1 | 0 | 4 | 0 | | | | |
| CD3D | 2 | 0 | 1 | 0 | | | | |
| CD3E | 2 | 0 | 2 | 0 | | | | |
| CD3G | 1 | 0 | 7 | 0 | | | | |
| CD8B | 3 | 0 | 3 | 0 | | | | |
| IL7R | 1 | 0 | 5 | 0 | | | | |
| LCK | 2 | 0 | 5 | 0 | | | | |
| LDHB | 1 | 0 | 2 | 0 | | | | |
| LEF1 | 1 | 0 | 8 | 0 | | | | |
| RPL10 | 1 | 0 | 9 | 0 | | | | |
| Previous | 1 | 2 | 3 | 4 | 5 | ... | 28 | Next |

This marker module will show Differentially Expressed genes in the following:

- 1. Each group (using findMarker in scan)
- 2. Each group and pathway-specific (using Fisher’s exact test)

Developers

- Hai Nguyen (hainct@unist.ac.kr)
Author, maintainer
- Jinhwan Kim
Author
- Dougu Nam (dougnam@unist.ac.kr)
Author, project manager

License

- [MIT](#) + file [LICENSE](#)